## APPENDIX 1 EQUIPMENT FOR FISTULA REPAIR

ESSENTIAL SPECIALIST EQUIPMENT FOR FISTULA REPAIRS		
Table	Table with steep head down tilt (Trendelenburg) of 45 degrees and where the foot end of the table can be removed as otherwise it will be very difficult for the surgeon to sit for any duration.	
	• If the foot end of the table is fixed and cannot be removed, check if you can turn the table around 180 degrees so that the head becomes the foot.	
	• If there is not adequate head-down tilt, you can place bricks or step under the end of the bed.	
	Tables of good quality are made by Fazzini, Nuova.	
Shoulder pads	These allow steep head down tilt without the patient slipping down the table (Fig.1.3c).	
Light source	A bright light source e.g. Evolution LED lights, Astramax	
_	• A headlight is useful, as the table position will be changing during the surgery.	
Suction	This is important to avoid using a lot of gauze for swabbing.	
Ureteric	Size 4, 5 and 6 F. If not available, infant feeding tubes can be used but are more difficult to	
catheters	keep in place.	
Sutures	Strong needle/ suture: for taking the lateral bites of para-urethral tissue including periosteum.	
needles The	A 5/8 circle 26 or 27 mm needle (available as UR6 Vicryl 2/0 J602H or Polysorb UL-878 or	
main sutures	Novosyn 2/0 FR26).	
for repairs:	<ul> <li><i>Fine needle/ suture:</i> when taking the central bites of the urethra: A half circle 26 mm needle Vicryl 2/0 taper (VCP or J 417).</li> <li><i>Very fine needle/ suture</i> e.g. work on a fragile urethra or ureter:</li> </ul>	
	• 3/0 22 mm needle (Polysorb GL-182) or Novosyn HR22 or Vicryl 3/0 J305 3/0 or	
	• Vicryl 2/0 VCP306 on a 17 mm 1/2c needle.	
	Special: A J shaped 22 mm needle (Vicryl W9223) is useful when working in confined spaces	
	as it is easier to turn the needle e.g. when the fistula angle is very lateral.	
	Vagina: any large needle e.g. 36 mm Vicryl J375H or Polysorb 40 mm CL-803.	
Cautery	Cautery (Bovie) is very useful especially for higher/ cervical fistulae.	



## Fistula repair set:

*Upper row:* towel clips, kidney dish with gallipots, sponge forceps.

*Middle row:* weighted speculum, Sims, Retractors, long and short artery forceps, long and short Allis forceps, vulsellum.

*Lower row:* Uterine sound, metal catheter, probe, dissecting forceps/ pickups, scissors, needle holders. On left, tray with suction tubing. Thorek scissors which are useful in fistula surgery are not shown in this set.

LIST OF INSTRUMENTS FOR FISTULA REPAIR SET			
Speculum + Retractors	<ul> <li>Auvard weighted speculum with an acute angle which allows it to sit nicely in the vagina. There are two sizes: (a) The larger one (Aesculap ® EL 739R) is suitable for most cases. (b) A smaller one with a shorter and narrower blade (Aesculap® EL 736R) is suitable for small or scarred vaginas and is particularly useful to have.</li> <li>Sims speculum: is useful for a shorter vagina.</li> <li>Deaver retractors:         <ul> <li>Small curved 8 ½ inch (Aesculap® BT601R) is very useful to help in retracting the bladder during large fistula repairs.</li> <li>Large retractors are useful in high fistulas especially in obese patients.</li> <li>Langenbeck retractor, 13 x 44 mm blade x 2</li> </ul> </li> </ul>		
Scissors	<ul> <li>Dissecting scissors, curved: Sharp ones are ideal e.g. Stille-Matarasso are particularly good for cutting scar tissue. The Stille tonsillar dissecting scissors are finer and sharp.</li> <li>Suture scissors, straight, 18 cm</li> <li>Thorek scissors</li> </ul>		
Needle holder	<ul> <li>Mayo-Hegar, 20.5 cm x 1</li> <li>Mayo-Hegar, 18 cm x 1</li> </ul>		
Dissecting forceps	<ul> <li>Dissecting forceps, with teeth, 20 cm x 1</li> <li>Dissecting forceps, fine serrated jaw, 20 cm x 1</li> </ul>		
Single items	<ul> <li>Blade holder 7, Swann Morton, 159 mm x 1</li> <li>Probe with eye, malleable, 20 cm x 1</li> <li>Uterine sound, malleable, 30 cm x 1</li> <li>Metal catheter, 16 cm x 1</li> <li>Kidney dish, large, 34 cm and gallipots x 2</li> </ul>		
Allis forceps	<ul> <li>Allis forceps, 20 cm x 2</li> <li>Allis forceps, 15 cm x 2</li> </ul>		
Artery forceps	<ul> <li>Large artery forceps, 23 cm x 1</li> <li>Mosquito's, curved, 13 cm x 4</li> </ul>		
Large forceps	<ul> <li>Sponge holding forceps, 241 mm x 2</li> <li>Vulsellum forceps, curved, 229 mm x 1</li> </ul>		
Miscellaneous	<ul> <li>Towel clips, Backhaus, 89 mm x 4</li> <li>Towel clips, Backhaus, 127 mm x 6</li> <li>Safety Pin forceps holder, 114 mm x 2</li> </ul>		
Special blades	Either size 11 or 15 are useful for incising in the vagina.		

CHECK LIST OF ITEMS TO USE FOR FISTULA REPAIRS		
Ureteric catheters		
Foley catheters sizes 12F 14F 16F		
Syringes: catheter tipped 60-100 ml		
Dye e.g. gentian violet paint or methyle	ne blue	
Surgical blades: size 11 or 15		
Special sutures		
Probes: to help catheterize ureters		
Prophylactic antibiotics		
Spinal anaesthesia e.g. Heavy Bupivaca	ine	
Post-operative instruction sheet (Appen	dix 3) and post-operative charts (Appendix 4).	
Strapping of good quality		
The 3 essential instruments: Good sciss	ors, needle holder and toothed dissecting forceps (pickups).	